


# Overview on data correction for polarized SANS

- (1) Overview  $^3\text{He}$ -cells and Supermirrors
- (2) Data reduction steps
  - (a)  $^3\text{He} + ^3\text{He}$  (DREAM)
    - DB measurements with 2 channels on both cells
  - (b) SM +  $^3\text{He}$  (SKADI)
    - DB measurements with 4 channels
    - DB measurements with 2 channels

# (i) Neutron polarizers and analyzers



(1) How are neutrons polarized?

→ Polarizer, "P"

(2) How are they transported?

→ Guide field, "G"

How are they flipped?


→ Flipper, "F"

How is the neutron polarization read out?

→ Analyzer, "A"

# (i) Neutron polarizers and analyzers: $^3\text{He}$ glass cells

Based on spin-dependent neutron absorption cross-section of  $^3\text{He}$




$^3\text{He}$ -cell:

- Filled with polarized  $^3\text{He}$  atoms
- Transmission of neutrons with  $\mathbf{s}_n \parallel \mathbf{s}_{\text{He}}$ 
  - $\sigma_{++} \approx 5$  barn
- Adsorption of neutrons with  $\mathbf{s}_n \parallel -\mathbf{s}_{\text{He}}$ 
  - $\sigma_{+-} \approx 6000$  barn
- Reversal of  $\mathbf{s}_{\text{He}}$  by NMR-pulses
- **Atomic polarization of  $^3\text{He}$ -atoms decays with time!**


$$P(^3\text{He})(t) = \tanh(\mu * AP(t))$$

$AP(t)$   
 = Atomic polarization  
 =  $AP_0 * \exp(-\Delta t / \gamma)$   
 $\mu$  = opacity of  $^3\text{He}$ -cell



$$T(^3\text{He})(t) = TE * \exp(-\mu) * \cosh(\mu * AP(t))$$

$TE$  = transmission of unfilled cell  
 $\gamma$  = decay time constant



Correction of measurements



$$I = I(\text{measured})(t) / \text{Trans}(t)$$



[1] K. Krycka et al., *Physica B* **404**, 2561-2564 (2009).  
 [2] R. P. Hermann, *Neutron Scattering Lectures*, Forschungszentrum Jülich GmbH (2015), vol. 106

# (i) Neutron polarizers and analyzers: Supermirrors




$\rho_{tot} = \rho_n \pm \rho_m$   
 $\rho(Fe, Spin+) = 12.97 \times 10^{-6} \text{ \AA}^{-2}$   
 $\rho(Fe, Spin-) = 3.08 \times 10^{-6} \text{ \AA}^{-2}$   
 $\rho(Si) = 2.08 \times 10^{-6} \text{ \AA}^{-2}$

Based on spin-dependent neutron scattering length of magnetized Fe layers!



### Supermirror:

- FeSi multilayers
- High reflectivity for one spin state, high transmission for other
- Broad total reflection niveau due to multiple Bragg angles
- „m-value“  $\approx$  number of layers  $\approx \theta_{tot,max}$



[1] S. Broekhuijsen. PhD thesis. University of Linköping (2021).

# (i) Neutron polarizers and analyzers: Supermirrors







Figure 3.7: While intermixing and surface roughness are distinct phenomena, they are indistinguishable in the SLD profile normal to the interface.

[1] S. Broekhuijsen, PhD thesis, University of Linköping (2021).

# (i) Neutron polarizers and analyzers: Supermirrors



$\rho_{tot} = \rho_n \pm \rho_m$   
 $\rho(Fe, Spin+) = 12.97 \times 10^{-6} \text{ \AA}^{-2}$   
 $\rho(Fe, Spin-) = 3.08 \times 10^{-6} \text{ \AA}^{-2}$   
 $\rho(Si) = 2.08 \times 10^{-6} \text{ \AA}^{-2}$



Different geometries:


Single mirror




V-cavity



# (ii) Example setup: Pol-SANS (SM + 3He)




# (ii) Example setup: PND ( $^3\text{He}$ + $^3\text{He}$ )






# (vi) Example setup: PNR (SM + SM)


A



B



# (iii) Data analysis under consideration of all components



Effect on measured intensities?

$$\begin{pmatrix} I_{++} \\ I_{+-} \\ I_{-+} \\ I_{--} \end{pmatrix} \sim \frac{d\sigma}{d\Omega} = \begin{pmatrix} S_{++} \\ S_{+-} \\ S_{-+} \\ S_{--} \end{pmatrix}$$

### (iii) Data analysis under consideration of all components

$$I \sim \frac{d\sigma}{d\Omega} = \begin{pmatrix} S_{++} \\ S_{+-} \\ S_{-+} \\ S_{--} \end{pmatrix} \quad \begin{pmatrix} I_{++} \\ I_{+-} \\ I_{-+} \\ I_{--} \end{pmatrix} = I_0 \cdot \begin{pmatrix} A_+ & A_- & & \\ A_- & A_+ & & \\ & & 0 & \\ & & A_+ & A_- \\ & & A_- & A_+ \end{pmatrix} \cdot \mathbf{F2} \cdot \mathbf{F1} \cdot \begin{pmatrix} P_+ & 0 & P_- & 0 \\ 0 & P_+ & 0 & P_- \\ P_- & 0 & P_+ & 0 \\ 0 & P_- & 0 & P_+ \end{pmatrix} \cdot \begin{pmatrix} S_{++} \\ S_{+-} \\ S_{-+} \\ S_{--} \end{pmatrix}$$

$I$  = measured intensities

$S$  = The sample's scattering cross-section  $\rightarrow$  what we want to get from the data reduction.

$P$  = Polariser spin-dependent transmission matrix

$A$  = Analyser spin-dependent transmission matrix

$\mathbf{F1}, \mathbf{F2}$  = spin-flipper matrices, if present  $\rightarrow$  If they work with 100% efficiency:  $\mathbf{F1}, \mathbf{F2} = \mathbf{1}$  (*identity matrix*)

$[P, A] = 0 \rightarrow P$  &  $A$  commute, order does not matter

$[P, F] \neq 0 \rightarrow$  the matrices must be in order!

# Background information on $^3\text{He}$ -cell transmission correction

## Definitions:

Nuclear polarization of  $^3\text{He}$  cell

$$P_{\text{He}}(t)$$

→ same for +/- states,  $\lambda$ -independent

Opacity of  $^3\text{He}$ -cell:

$$O(\lambda)$$

→ same for +/- states,  $\lambda$ -dependent

Transmission empty  $^3\text{He}$  glass:

$$TE$$

→ same for +/- states,  $\lambda$ -independent

Transmission unpolarized neutron beam:

$$Transmission^{unpol}(t, \lambda)$$

→ same for +/- states,  $\lambda$ -dependent

## Corrections to be done for each He cell:

### 1. Opacity + unpolarized cell transmission

(a) Either by measurement of pressure from  $^3\text{He}$ -team after filling:

$$O(\lambda) = 0.0733 p l \lambda$$

(b) Or by neutron beam:

$$O(\lambda) = -\ln \left[ \frac{I_{unpolarized\ ^3He\ cell}^{unpol.beam} - I_{bg, He}}{I_{^3He\ cell\ out}^{unpol.beam} - I_{bg, noHe}} \frac{1}{TE} \right]$$

### 2. Initial nuclear polarization

$$P_{\text{He}0}(t_0) = \text{acosh} \left[ \frac{I_{polarized\ ^3He\ cell}^{unpol.beam} - I_{bg, He}}{I_{^3He\ cell\ out}^{unpol.beam} - I_{bg, noHe}} \frac{1}{TE \cdot \exp(O(\lambda))} \right] / O(\lambda)$$

### 3. Time-decay

$$P_{\text{He}}(t) = P_{\text{He}0}(t_0) * \exp\left(-\frac{\Delta t}{T_1}\right)$$

### 4. Determine correction factors for polarization dependent cell transmission

$$Transmission^{\pm}(t, \lambda) = TE * \exp(-O(\lambda) \pm O(\lambda) * P_{\text{He}}(t))$$

Details are listed for each step on the next slides

# Detailed Data-correction procedure

## 1. Opacity + unpolarized cell transmission

There are 2 ways to measure the opacity:

(a) By measurement of pressure from  $^3\text{He}$ -team after filling

$$O(\lambda) = 0.0733 p l \lambda$$

- disadvantage: measurement of p not very precise
- advantage: Can work as a fast first estimation without measurement of the unpolarized cell

<b>Fix input values</b>	Cell pressure, p	Cell length, l	Wavelength, $\lambda$
<b>Neutron data input</b>	None		
<b>Parameters to compute</b>	Opacity O ( $\lambda$ )		

(b) By neutron beam measurements:

$$O(\lambda) = -\ln\left[\frac{I_{\text{unpolarized } ^3\text{He cell}}^{\text{unpol}} - I_{\text{bg, He}}}{I_{^3\text{He cell out}}^{\text{unpol}} - I_{\text{bg, noHe}}} \frac{1}{TE}\right]$$

- advantage: more accurate
- disadvantage: cell needs to be depolarized, i.e., would be done only at the end of its lifetime (i.e., after ca. 1 day)
- (b) Should be done at the end of the day/experiment for very precise data reduction
- (a) should be taken in the beginning for live-data reduction, but afterwards being overwritten by (b)

<b>Fix input values</b>	Transmission empty glass, TE	Wavelength, $\lambda$
<b>Neutron data input</b>	DB transmission of - unpolarized neutron beam - unpolarized $^3\text{He}$ cell $I_{\text{unpolarized } ^3\text{He cell}}^{\text{unpol}}$	DB Transmission of - unpolarized neutron beam - no cell in (moved out) $I_{^3\text{He cell out}}^{\text{unpol}}$
<b>Parameters to compute</b>	Opacity O ( $\lambda$ )	

# Detailed Data-correction procedure

## 2. Initial nuclear polarization

**Step Nr. 2 : Determination of time-independent prefactor of the  $^3\text{He}$ -cell time decay:**

$$P_{He0} = \text{acosh} \left[ \frac{I_{\text{polarized } ^3\text{He cell}}^{+ \text{ or } -} - I_{bg, He}}{I_{^3\text{He cell out}}^{unpol} - I_{bg, noHe}} \frac{1}{TE \cdot \exp(O(\lambda))} \right] / O(\lambda)$$

**To note:**

**(1)**  $P_{He0}(t_0)$  is polarization independent

- Up/ Down polarized cells would give same result for the direct beam measurement:  $I_{\text{polarized } ^3\text{He cell}}^{+} = I_{\text{polarized } ^3\text{He cell}}^{-}$

- Take only one polarization direction for the calculations (saving time)

**(2)**  $T_{bg}$  is to be taken from the same measurement as  $T$ , but far away from Direct Beam region

- probably need some maskfile

- this is important in case the background with  $^3\text{He}$  cell is different than without  $^3\text{He}$  cell

<b>Fix input values</b>	Transmission empty glass, TE	Opacity, $O(\lambda)$
<b>Neutron data input</b>	DB transmission of - half-polarized neutron beam - polarized $^3\text{He}$ cell (any pol.-direction) - according background data $I^{+}$ or $I^{-}$	DB Transmission of - unpolarized neutron beam - no cell in (moved out) - according background data $I_{^3\text{He cell out}}^{unpol}$
<b>Parameters to compute</b>	Nuclear polarization at time $t_0$ , $P_{He0}(t_0)$	

# Detailed Data-correction procedure

## 3. Time-decay of nuclear polarization:

**Step Nr. 3 : Determination of time-dependency of each  $^3\text{He}$ -cell:**

$$P_{\text{He}}(t_n) = P_{\text{He}0} * \exp\left(-\frac{t_0 - t_n}{T_1}\right)$$

**To note:**

**(1)**  $P_{\text{He}}(t)$  is polarization independent

- Up/ Down polarized cells would give same result

- Take only one polarization direction for the calculations (saving time)

**(2)**  $T_{\text{bg}}$  is to be taken from the same measurement as  $T$ , but far away from Direct Beam region

- probably need some maskfile

- this is important in case the background with  $^3\text{He}$  cell is different than without  $^3\text{He}$  cell

<b>Fix input values</b>	Start time of using $^3\text{He}$ cell in the beamline (start of first measurement), $t_0$	Time stamp $t_n$
<b>Neutron data input</b>	DB transmission of - unpolarized neutron beam - polarized $^3\text{He}$ cell (P+A) (any pol.-direction) - take this as function of time!  $I^+$ or $I^-$	
<b>Parameters to compute</b>	Nuclear polarization at time $t_n$ , $P_{\text{He}}(t_n)$	Decay time $T_1$

# Background information on $^3\text{He}$ -cell transmission correction

## 4. Calculation of correction parameters for the sample measurements:

**Step Nr. 4 :** is to define the polarization dependent transmission functions of each  $^3\text{He}$  cell (Polarizer, P, and analyzer, A):

$$T^\pm(t_n, \lambda) = \text{TE} * \exp(-O(\lambda) \pm O(\lambda) * P_{\text{He}}(t))$$

$$\begin{pmatrix} I_{++} \\ I_{+-} \\ I_{-+} \\ I_{--} \end{pmatrix} = I_0 \cdot \begin{pmatrix} A_+ & A_- & & \\ A_- & A_+ & & \\ & & A_+ & A_- \\ & & A_- & A_+ \end{pmatrix} \cdot \begin{pmatrix} P_+ & 0 & P_- & 0 \\ 0 & P_+ & 0 & P_- \\ P_- & 0 & P_+ & 0 \\ 0 & P_- & 0 & P_+ \end{pmatrix} \cdot \begin{pmatrix} S_{++} \\ S_{+-} \\ S_{-+} \\ S_{--} \end{pmatrix}$$

**To note:**

**(1) The transmission is only polarization dependent for an incoming polarized neutron beam, hence polarization dependence of the transmission needs to be taken into account only for the sample measurements, but not for determining the time-decay of each single  $^3\text{He}$ -cell**

Fix input values	Transmission of empty glass, TE	Opacity, $O(\lambda)$		Time-stamp, $t_n$		Wavelength, $\lambda$	
Neutron data input	Sample data of - full-polarized neutron beam, all four spin-channels - $I^{++}, I^{+-}, I^{-+}, I^{--}$						
Parameters to compute	Polarization-dependent transmissions through the cell - $T_+, T_-$ (for polarizer and analyzer)	Corrected sample data - $I^{++}_{\text{corr}}, I^{+-}_{\text{corr}}, I^{-+}_{\text{corr}}, I^{--}_{\text{corr}}$					



# Important information

---

1. Keep Wavelength-dependence for PA-correction
2. Normalization
  - normalization to monitor counts should be done before PA-correction
  - probably good to also normalize to all other background noise / detector efficiencies, etc. before PA-correction
  - think about: position of Monitors (for analyzer before, but for polarizer after cell?) → does that make a difference?
3. Difference of SM /  $^3\text{He}$ -cell usage:
  - Use same reduction script for any  $^3\text{He}$ -cell, independent of choice of polarizers
4. To think about:
  - 3D detector for SANS at DREAM – ask Celine/Wojtek/Mikhail how this would work
  - DB background same for using the  $^3\text{He}$ -cell or not? Is DB background polarization dependent?
  - Have to test this prior to operation, but consider it as same for now

# Workflow for $^3\text{He} + ^3\text{He}$ – Option 2: 2 channels (W/o polarizer) + unpol

Measure parameters of both  $^3\text{He}$  cells

→ Opacity, pressure, initial polarization,...

→ To be fitted afterwards:



Sample: Measure 4-channels with sample

→ (w polarizer, w sample, w analyzer)

→ Do 15 min measurements (in our case)



$^3\text{He}$ :

(i) Measure unpolarized beam

→ Everything out, reference measurement

(ii) Measure I+ or I- (DB) ( $^3\text{He}1$ )

→ (without polarizer)

(iii) Measure I+ or I- (DB) ( $^3\text{He}2$ )

→ (without analyzer)

(iv) as 60s measurements

→ Good time-steps: 1h

E.g. following procedure

- 4-channel measurement on sample
  - Move sample in
  - Measure 15 mins x 4 channels (ca 1h loop)
- Unpol-measurement on DB ("unpol")
- I+ or I- measurement on DB ("pol out")
  - Move polarizer out
  - Measure 60 s x 2 channels
- I+ or I- measurement on DB ("an out")
  - Move analyzer out
  - Move polarizer in
  - Measure 60 s x 2 channels

→ **Loop** until good signal to noise ratio for the measurements

For future: workflows for Supermirror +  $^3\text{He}$  cell

# Workflow for SM + 3He – Option 1: 4 channels – **to be optimized**

$$\begin{pmatrix} I^{++} \\ I^{+-} \\ I^{-+} \\ I^{--} \end{pmatrix} = I_0 \begin{pmatrix} A^+ & A^- & 0 & 0 \\ A^- & A^+ & 0 & 0 \\ 0 & 0 & A^+ & A^- \\ 0 & 0 & A^- & A^+ \end{pmatrix} \begin{pmatrix} 1 & 0 & 0 & 0 \\ 0 & 1 & 0 & 0 \\ 1 - F1 & 0 & F1 & 0 \\ 0 & 1 - F1 & 0 & F1 \end{pmatrix} \begin{pmatrix} P^+ & 0 & P^- & 0 \\ 0 & P^+ & 0 & P^- \\ P^- & 0 & P^+ & 0 \\ 0 & P^- & 0 & P^+ \end{pmatrix} \begin{pmatrix} S^{++} \\ S^{+-} \\ S^{-+} \\ S^{--} \end{pmatrix}$$

$P(t, \lambda)$  = Polarizer transmission

$A(t, \lambda)$  = Analyzer transmission

$F(\lambda)$  = Flipper efficiency

(1) Normalize data by counting time or the monitor counts.

(2) Functions:

$$A^\pm(t, \lambda) = TE * \exp(-O(\lambda) \pm O(\lambda) * P_{He}(t))$$

$A(t, \lambda)$  = Analyzer (3He) transmission

$$P_{He}(t) = P_{He0}(t_0) * \exp\left(-\frac{\Delta t}{\gamma}\right)$$

AP(t) = Atomic polarization

## Known Values:

Opacity	$O(\lambda) = 0.0733 p l \lambda$
filling time	$t_0$
initial He-polarization	$P_{He0}(t_0)$
time-stamp	$t_n$

## Unknown params to fit:

Base Intensity	$I_0$
Flipper efficiency:	$E_{F1}(\lambda)$
3He time-decay constant:	$T1$ [h]
Transmission unfilled cell:	TE

$P$   
 $\lambda$   
 $l$   
pressure of 3He cell  
neutron wavelength.  
path length through the cell

## Can be either pre-measured or fitted:

Polarizer transmission SM(+):	$P^+(\lambda)$
Polarizer transmission SM (-):	$P^-(\lambda) = 1 - P^+(\lambda)$

Measure  $I$   $\rightarrow$  Fit  $A^\pm(t, \lambda)$   $\rightarrow$  Insert into Matrix

# Workflow for SM + $^3\text{He}$ (MEOP) – Option 1: 4 channels

Measure parameters of  $^3\text{He}$  cell

→ Opacity, pressure, initial polarization,...

→ To be fitted afterwards:



Sample: Measure 4-channels with sample

→ (w polarizer, w sample, w analyzer)

→ Do 15 min measurements (in our case)



$^3\text{He}$ : Measure 4 channels without sample

→ (w polarizer, w analyzer)

→ Do 60s measurements

→ Max. time-steps: 2h (what we did at D33)

→ Good time-steps: 1h

E.g. following procedure (from D33 May 2023 Annika):

- 4-channel measurement on sample
  - Move sample in
  - Measure 15 mins x 4 channels (ca 1h loop)
- 4-channel measurement on DB
  - Move sample out
  - Measure 60 s x 4channels

→ **Loop** until good signal to noise ratio for the measurements

## Additionally to consider:

- Sample depolarization – effectively changes polarizer efficiency, but would still be nice to be able to insert if users know it?
  - Proposition: leave it out now, it would just come in as another matrix element with one additional factor
- At how many wavelength to measure?
  - Binning DB = sample binning

# Workflow for SM + $^3\text{He}$ – Option 2: 2 channels + unpol – **to be optimized**

## Input Values:

filling time	$t_0$
time-stamp	$t_n$
Wavelength	$\lambda$
Polarizer transmission SM(+):	$P^+(\lambda)$
Polarizer transmission SM (-):	$P^-(\lambda) = 1 - P^+(\lambda)$

## Unknown params to fit:

Base Intensity	$I_0$
$^3\text{He}$ time-decay constant:	$\gamma$ [h]
Transmission unfilled cell:	TE

(a) Either by measurement of pressure from  $^3\text{He}$ -team after filling:

$$O(\lambda) = 0.0733 p l \lambda$$

(b) Or by neutron beam:

$$O(\lambda) = -\ln \left[ \frac{T_{\text{unpolarized } ^3\text{He cell}}^{\text{unpol.beam}} - T_{\text{bg, He}}}{T_{^3\text{He cell out}}^{\text{unpol.beam}} - T_{\text{bg, noHe}}} \frac{1}{TE} \right]$$

## 2. Initial nuclear polarization

$$P_{\text{He0}}(t_0) = \text{acosh} \left[ \frac{T_{\text{polarized } ^3\text{He cell}}^{\text{unpol.beam}} - T_{\text{bg, He}}}{T_{^3\text{He cell out}}^{\text{unpol.beam}} - T_{\text{bg, noHe}}} \frac{1}{TE \cdot \exp(O(\lambda))} \right] / O(\lambda)$$

## 3. Time-decay

$$P_{\text{He}}(t) = P_{\text{He0}}(t_0) * \exp\left(-\frac{\Delta t}{T_1}\right)$$

## 4. Determine correction factors for polarization dependent cell transmission

$$\text{Transmission}^\pm(t, \lambda) = TE * \exp(-O(\lambda) \pm O(\lambda) * P_{\text{He}}(t))$$

# Workflow for SM + $^3\text{He}$ – Option 2: 2 channels (W/o polarizer) + unpol

Measure parameters of  $^3\text{He}$  cell

→ Opacity, pressure, initial polarization,...

→ To be fitted afterwards:



Sample: Measure 4-channels with sample

→ (w polarizer, w sample, w analyzer)

→ Do 15 min measurements (in our case)



$^3\text{He}$ :

(i) Measure unpolarized beam

→ W/o polarizer, w/o analyzer, w/o sample

(ii) Measure  $I^+$  or  $I^-$

→ (without polarizer, with analyzer)

(iii) Loop as 60s measurements

→ Good time-steps: 1h

E.g. following procedure

- 4-channel measurement on sample
  - Move sample in
  - Measure 15 mins x 4 channels (ca 1h loop)
- Unpol measurement on DB (Tunpol,  $^3\text{He}$  out)
  - Move polarizer out
  - Move analyzer out
  - Measure 60 s
- $I^+$  or  $I^-$  measurement on DB ( $^3\text{He}$  in)
  - Move analyzer in
  - Measure 60 s x 2 channels

→ **Loop** until good signal to noise ratio for the measurements